



**HAEBERLE**

mobile workstations

**HAEBERLE  
MAKES IT**

# FIRST-CLASS TRANSPORT FOR YOUR EQUIPMENT.


Multifunctional cart systems with added value.

HAEBERLE sounds Swabian and is Swabian: inventive, technophilic, diligent, reliable; in short, a maker of unique customised carts. We manufacture equipment cart systems for many sectors, but first and foremost for medical technology. Our systems are made in Stuttgart and renowned all over the world.

What makes us special? Our carts offer first-class functionality and aesthetics. After all, not even the best equipment can be used efficiently if the cart system lets it down. With us, your equipment can be transported on a cart with all the extras.

Perhaps one with pedal mounts or docking points for additional equipment, laptop and monitor brackets and much more besides – there are virtually no limits. The same goes for the design: we design the cart systems to match your corporate identity and enable your equipment to be perfectly integrated.


Since 1933, we have been coming up with the perfect solution for each and every customer; solutions in line with all of the technical rules and standards, in an uncomplicated and reliable fashion. **Sound good to you? Well, it is! Find out for yourself; we look forward to hearing your special wishes.**



# 1 BASIC

## Components, accessories, colours, logo placement

A brand should be instantly recognisable, which is why we adapt our cart systems to mirror your corporate identity. The important aspects here are the fittings, colour coordination and placement of your logo in the right position.




**MONITOR HOLDER** „alegro“ with fixing bracket for stepless mounting to the carrier profile.



**LIFT-FUNCTION** for sitting or standing working position by adjusting the smooth movable carrier unit.



**TRIPPING LEVER** for stepless mechanical height adjustment.



**PLATES** made of ABS plastic, steel plate or here for example with a keyboard and mouse pull-out shelf.



**ACCESSORIES**, such as here, a PC support for fastening to the carrier profile or other accessory parts are possible.



**COLOURS** to choose for decorative strips + logo of your company, for example, on the chassis.

### **LITTLE EFFORT, GREAT EFFECT:**

Varying the basic components alone allows you to individually design your cart.

For example the “fuego” cart with comfortable lift function

# 2

## BASIC PLUS

### Basic models with special upgrades

Each person is different, so too is each working area. This is why a system is only good if it suits the individual application. Whether you have plenty of space or work in very tight confines, whether you need a laptop cart or an ECG cart or something else entirely, we are certain to find the solution.

#### ENDLESS POSSIBILITIES FOR ALL FUNCTIONS:


- Dimensions
- Colours
- Special shelves
- Pull-out shelves
- Computer
- Keyboard arms
- Endoscopy
- Ergometry
- Chassis
- Articulated arms
- Handles
- Cable channel
- Monitor holder
- Infusion
- Emergency
- Double monitoring
- Carrier profile
- Castors
- Drawers
- Sonography
- Stroboscopy
- Lift functions
- custom imprinting options
- ... and more



"fuego Plus" for the CURIS RF generator from Sutter Medizintechnik GmbH with holding arm for the pedal and baskets for accessories.



"bravo Plus" Dental cart has been designed in cooperation with KaVo Dental GmbH specifically for the surgical units MASTERSurg, EXPERTSurg and INTRASurg.



"fuego Plus" the special cart for the Touch device of LASERNEEDLE GmbH for orthopedics and sports therapy.

## A CASE IN POINT

**THE CHALLENGE:** The Italian company EBneuro desired a functional cart for its "GastroExplorer" equipment for examination of the oesophagus.

**PRECONDITION:** Water and power connections were required to facilitate the diverse functions. In addition, the system had to be equipped with a connection for a computer with an easily visible monitor. To protect all the equipment, stability was a crucial factor in component selection.

**SOLUTION:** A toro cart specially tailored for the equipment.



imt-medical ag, Switzerland, employs a fuego Plus for bellavista, its mobile display for patient data, ventilator settings and much more.

# 3 DEVELOPMENT

New outline to your specifications

Do you need a special cart for a specific application – tailor-made, if you will? You've come to the right place: From idea to consultation, design, prototype and on to production, we develop the customised cart with you.

## Step by step to your customised cart

### \_ Development consultation

Consulting relies on trust. Whatever your concerns about the development and design of your cart system, we are at your disposal.

1

### \_ Order

After the consultation, you decide from the options available to you and place an order with the specifications you require.

2

### \_ Designs, visualisations

We present you with concrete designs and 3D visualisations.

3

### \_ Models

As a check for the final version, we create models to test functionality, ergonomics, technology and design.

4

### \_ Construction


Once the drafts have been accepted and the model approved, it is the job of the engineering department to skilfully implement the specifications.

5

### \_ Prototypes

We create prototypes based on the design documents to get the go-ahead for production.

6





# HAEBERLE

mobile workstations

**HAEBERLE  
HAS IT**

**HAEBERLE  
MAKES IT**

HAEBERLE HAS IT: a modular basic system for your multifunctional carts. And if you need something extra? No problem, remember HAEBERLE MAKES IT: exactly the customised solution you require.

**HAEBERLE GmbH + Co. KG**

Breitwiesenstr. 13, DE-70565 Stuttgart  
PO Box 800524, DE-70505 Stuttgart

Phone +49 (0) 711 78314-0  
Fax +49 (0) 711 78314-77

[info@haeberle-med.de](mailto:info@haeberle-med.de)  
[www.haerberle-med.de](http://www.haerberle-med.de)

