

Dental Cart flexible and mobile in the OR

- Optimum compatibility, e.g. with the KaVo MASTER-/ EXPERTsurg & INTRAsurg surgical devices
- Optimum space utilisation due to large storage space on a small surface area
- Maximum mobility due to light weight and high-quality easy-run castors
- Ideal working position for activities whilst standing and sitting



Features: Dental Cart from HAEBERLE, Stuttgart

- The cart system stands out due to its modern design and stability and has been developed jointly with KaVo Dental GmbH
- Different storage places for the individual KaVo surgical devices components (see illustration)
- Smooth surfaces allow for optimum cleaning options
- High-quality shelf with shaped edge for ensuring the devices are secure when the cart is being manoeuvred
- "Made in Germany" – Guaranteed high quality
- Dental Cart is available in a choice of two colours (white/grey)



Chassis / white



Storage bowl / white



ABS shelf / grey

HAEBERLE **Equipment cart**

Item	Description	Mat. No.
Dental Cart white	Functional cart for KaVo surgical devices, decorative colour grey, ABS shelf on top, storage bowl white with frame and lid, foot starter bracket, chassis white	27097
Dental Cart grey	Functional cart for KaVo surgical devices, decorative colour grey, ABS shelf on top, storage bowl white with frame and lid, foot starter bracket, chassis grey	27099

Materials / choice of colour:

Carrier profile:	Anodised aluminium
Decorative strips in the profile:	Made of PVC in the grey colour version (6502B acc. to NCS)
Chassis:	Cast aluminium, powder-coated, grey (RAL 7012), fine structure, matt
Castors:	Polyurethane
Storage shelf (ABS):	ABS orion grey (surface: erosion structure 36)
Storage bowl:	Polycarbonate, white
Holder for storage bowl:	Steel, powder-coated (RAL 9006)
Cover for storage bowl:	Steel coated (RAL 9006)
Foot starter bracket:	Steel, powder-coated (RAL 9006)

Technical details / Dimensions:

Maximum load	Chassis:	25 kg	
	Storage shelf on top (ABS):	15 kg	
	Bowl:	5 kg	
	Foot starter bracket:	5 kg	
Dimensions	Chassis:	577 x 577 x 203 mm (W x D x H)	
	Carrier profile:	780 mm	
	Castors Ø:	75 mm, 2 lockable	
	Storage shelf on top (ABS):	External dimensions	500 x 390 x 53 mm (W x D x H)
		Storage surface	470 x 340 mm (W x D)
	Bowl:	External dimensions	280 x 210 x 80 mm (W x D x H)
		Internal dimensions	230 x 160 x 70 mm (W x D x H)
	Foot starter bracket:	External dimensions	240 x 280 x 78 mm (W x D x H)
Internal dimensions		189/229 x 272 mm (W1/W2 x D)	
Total weight	Dental Cart	13,5 kg	

Your dental distributor

KaVo is a brand name of the Kaltenbach and Voigt GmbH group of companies and is therefore a registered trade mark.