

HAEBERLE

mobile workstations

new in metal version

flexible and versatile

fuego equipment cart system
light duty



fuego
by HAEBERLE

fuego

The flexible, versatile equipment cart system

Not just for small spaces:
there's plentiful storage space over a small area

In the medical sector, fuego equipment carts are ideally suited for accommodating small equipment combinations in a mobile and flexible manner. They are also suitable for use in all types of offices and meeting rooms where the elegant fuego carts offer mobile solutions. In addition to high functionality, the exceptional design and extraordinary colours of fuego equipment carts are also captivating. They are versatile but also very sturdy and offer plentiful storage with a small footprint.

Customised is our standard – We consistently implement the basic modular idea to develop practical solutions. The wide variety of standard components makes each individual cart an individually tailored customer solution – whether it's for a doctor's practice or small or large quantities for equipment manufacturers. Even if you need to have it tailored just for you – HAEBERLE will make this possible.



fuego is protected with a European design patent and is a symbol for reliable brand quality at a reasonable price.



The system components

- Sturdy and very versatile chassis frame with high-grade, quality castors. The front castors are lockable and fitted with wall guards.
- Stable single-support aluminium profile for multiple socket outlets, with large integrated cable channel.
- High-quality trays made of plastic with shaped edge and recessed grip or metal.
- Special sizes available on request.
- Shapely drawers made of plastic or metal; easy-glide spherical telescopic guide; several drawers can be inserted under a plate.
- Produced as standard with large, ergonomic push handle; high-grade powder coating in RAL 9006 white aluminium.
- 3 decorative colours; also available in special colours for an extra charge.
- Accessories for many fields of application.

Equipment carts

Basic carts



fuego 1

Ref. no. 3000.1



fuego 1 with drawer

Ref. no. 3001.1



fuego 2

Ref. no. 3003.1

Computer carts



fuego - lift cart for notebook

Ref. no. 1978.1



fuego - computer cart 1

Ref. no. 1975.1



fuego - PC cart with lift

Ref. no. 1979.1

Ergometry cart / ECG carts



fuego - ergometry cart 1

Ref. no. 1962.1



fuego - ECG cart 2.1

Ref. no. 1961.1



fuego - ECG cart 1.1

Ref. no. 1960.1

Endoscopy cart / Sonography cart / Infusion cart



fuego - endoscopy cart

Ref. no. 3008.1



fuego - sonography cart

Ref. no. 3007.1



fuego - infusion cart

Ref. no. 3019.1

The fuego equipment cart system



- **new** with metal drawers and plates

Basic carts

new



fuego-metal 1

Ref.-Nr. 3024.1



fuego-metal 1 with drawer

Ref.-Nr. 3022.1



fuego-metal 2

Ref.-Nr. 3023.1

System parts metal

new



Metal drawer with metal plate

Ref. no. 3025.1



Undermount metal drawer

Ref. no. 3026.1



Steel plate

Ref. no. 19330

Choose from 3 colours for decorative strips and drawer handles.

In your order please specify the colour code by adding it as the last digit to the reference number:

6502B accord. to NCS



grey Ref. no. xxxx.1

3050-R80B accord. to NCS



blue Ref. no. xxxx.2

0050-Y accord. to NCS



yellow Ref. no. xxxx.4

Other colour combinations or individual colours on request.

System parts and accessories

HAEBERLE

mobile workstations



Drawer (ABS plastic)

Ref. no. 19300



Plate (ABS plastic)

Ref. no. 19033

new



Standard profile rail

Ref. no. 30016



Keyboard and mouse pullout

Ref. no. 19220



Laptop plate

Ref. no. 19320



Front handle for ABS plate

Ref. no. 19367



Filing basket in front

Ref. no. 19135



Keyboard and mouse pull-out shelf

Ref. no. 19222



"alegro" monitor holder

Ref. no. 19497



Infusion bottle holder

Ref. no. 19403



Cable and hose holder

Ref. no. 19034



Multiple standard socket outlet

Ref. no. 21696

Additional system components and accessories can be found in the fuego price list.

Your medical distributor

HAEBERLE GmbH + Co. KG

Breitwiesenstr. 13, DE-70565 Stuttgart

P.O.Box 800524, DE-70505 Stuttgart

Phone +49 (0) 711 78314-0

Fax +49 (0) 711 78314-77

info@haeberle-med.de

www.haerberle-med.de