

HAEBERLE

mobile workstations

versatile and good value

08/16[®] multipurpose trolley system



08/16[®]

08/16[®]

Tried-and-tested multipurpose trolley system

Just as the wide range of combinations impresses,
so too the versatility of applications

Our 08/16® white basic trolleys with colourful accents have been demonstrating their value in practice for some decades now and their infinite combination possibilities have won over users in a wide range of different application areas.

The trolleys are distinguished by their sturdy chassis with a white square-tubular, plastic-coated steel frame and are fitted with easy-care shelves and drawers with gentle-moving ball-bearing telescopic guidance and a click-in mechanism.

The trolley system is available as a simple basic trolley, an equipment tower, computer trolley, standing workstation, cabinet trolley ... and more. 08 /16® trolleys are used in doctors' practices, beauty salons, clinics and rehabilitation centres.

08/16®

08/16® – the multipurpose trolley system
from HAEBERLE – specialists in the
development of innovative cart systems
for medical equipment since 1933.



The system components

- 20 × 20 mm white square-tubular, plastic-coated steel frame.
- Gentle-running castors in 3 sizes, available on request with brakes or electrically conductive.
- Shelves with or without a raised edge.
- White light-resistant drawers and shelves, plastic-varnished on all sides.
- High-quality gently hovering ball-bearing telescopic drawer guidance with click-in mechanism.
- Ergonomic drawer handles available in 7 colours.
- 4 widths of trolleys 34 / 49 / 60 / 90 cm.
5 heights of trolleys 52 / 70.5 / 73 / 80.5 / 83 cm.
- Accessories for a wide range of applications.

Examples of model versatility

Basic trolleys



08/16® basic trolley 30
705 mm high

Ref. no. 15871



08/16® basic trolley 45
805 mm high

Ref. no. 15079



08/16® basic trolley 60
730 mm high

Ref. no. 12077



08/16® underneath trolley 60
520 mm high

Ref. no. 12180



08/16® basic trolley 90
830 mm high


Ref. no. 17277



08/16® cabinet trolley 60
830 mm high

Ref. no. 15862

08/16® basic trolleys



	width 34 cm	width 49 cm	width 60 cm	width 60 cm	width 60 cm	width 90 cm
number of drawers	height 70,5 cm	height 70,5 cm	height 52 cm	height 73 cm	height 83 cm	height 83 cm
1	15871	15139	12078	12075	14754	17277
2	10410	15155	12179	12076	14755	11008
3	10411	15156	12180	12077	14756	11020
4	10412	15041		17205	14757	
5					15391	

Shelf trolleys



08/16® shelf trolley 60 - mini
730 mm high / 260 mm deep

Ref. no. 12057



08/16® shelf trolley 60
830 mm high

Ref. no. 15150



08/16® shelf trolley 90
830 mm high

Ref. no. 17153

Function trolleys



08/16® standing working place 60
1000 mm high

Ref. no. 17120



08/16® ergometry trolley 60
910 mm high with console



Ref. no. 15167



08/16® injection trolley 60
1290 mm high with PicBox® Plus

Ref. no. 15154

08/16® self trolleys



	width 34 cm depth 45 cm	width 49 cm depth 45 cm	width 60 cm depth 26 cm	width 60 cm depth 45 cm	width 90 cm depth 45 cm
height 70,5 cm	15894				
73,0 cm			12057		
80,5 cm		14514			
83,0 cm				15150	17153

Customised is our standard

For maximum mobility in the workplace

Since 1933, our company has been employing smart ideas and applying practical solutions to develop countless innovative product lines in the area of multipurpose equipment carts. These have been marketed worldwide with great success.

With our products and solutions, we guarantee maximum mobility, flexibility and reliability in the workplace. This means everyone can focus fully and completely on their own work.








We consistently implement the basic modular idea to develop practical solutions. The wide variety of standard components makes each individual cart an individually tailored customer solution – whether it's for a doctor's practice, small volume production or large quantities for equipment manufacturers. Even if you need to have it tailored just for you – HAEBERLE will make this possible.



HAEBERLE HAS IT: a modular basic system for your multifunctional carts. And if you need something extra? No problem, remember HAEBERLE MAKES IT: exactly the customised solution you require.

Variety of colours

Choose from 7 fashionable colours for decorative handles

				
grey	blue	red violet	yellow	turquoise
Ref. no. 13161	Ref. no. 13162	Ref. no. 13163	Ref. no. 13164	Ref. no. 13165
				
white	black			
Ref. no. 13166	Ref. no. 13167			

Other colour combinations or individual colours on request.

System parts and accessories

HAEBERLE

mobile workstations



Drawer inset for trolley 60

Ref. no. 16723



Drawer rear panel for trolley 60

Ref. no. 15752



Waste collector „tomo“

Ref. no. 11026



Waste container

Ref. no. 04615



Injection set PicBox®

Ref. no. 11381



Injection set PicBox® multi

Ref. no. 16335



„Medi-Müll set“ 20 boxes 1,5 ltrs.

Ref. no. 04778



Ampoule opener „CUPFIX“

Ref. no. 05771



Arm support

Ref. no. 09496



Cable and hose holder

Ref. no. 10265



Infusion bottle holder

Ref. no. 04621



Console for trolley 60

Ref. no. 09499

Additional system components and accessories can be found in the 08/16® price list.

Your medical distributor

HAEBERLE GmbH + Co. KG

Breitwiesenstr. 13, DE-70565 Stuttgart

Postfach 800524, DE-70505 Stuttgart

Telefon +49 (0) 711 78314-0

Telefax +49 (0) 711 78314-77

info@haeberle-med.de

www.haerberle-med.de