

HAEBERLE

mobile workstations

variable and classy

keo® ISO module and
function trolley system



keo®

keo®

Variable and extraordinary trolley system

Endless scope for the creation of individual trolley system combinations

keo® is a modern and mobile trolley series with infinite version possibilities and clean lines with colourful accents, which is ergonomic, perfectly suited to its function and practical for many treatment applications. Select a suitable colour and height for your keo® function trolley, either with a metal work surface or with a raised-edge plastic plate, with premium-quality closable roller blinds or alternatively with drawers or ISO module full pull-outs, file pull-outs ... and more.

Various colour combinations are available with a distortion proof case made from high-quality, coated steel sheet. All sorts of accessories can be mounted on the four elegant aluminium support profiles.

The trolleys are suitable for use as dressing and treatment trolley, ISO module trolley, anaesthesia trolley and emergency trolley especially in clinics, rehabilitation centres and doctors' practices.

keo®

keo® – the ISO module and function trolley system from HAEBERLE – specialists in the development of innovative cart systems for medical equipment since 1933.



The system components

- Rugged, distortion proof metal construction made from galvanised, high-quality epoxy-resin steel sheet with a powder coating.
- 4 elegant anodised aluminium support profiles for mounting accessories. Subsequent mounting of accessories is also possible.
- All trolleys have high-quality smoothly running castors as standard. Ø 125 mm, with two brakes and with wall guards.
- All trolleys are fitted with an elegant push handle (white aluminium RAL 9006) as standard.
- Work surfaces optionally grey ABS plastic with raised edge, made of stainless steel with surge edge or powder-coated metal with lateral pull-out.
- All module full pull-outs have an inset frame for mounting conventional ISO modules. The pull-out can be extended up to approx. 110 %.
- Drawers can likewise be extended fully (110%).
- Decorative strips and drawer handles are available as standard in 5 colours.
- Special models available on request.

Examples of function trolleys

ISO module trolleys



ISO module trolley 1

Ref. no. 1851.1



ISO module trolley 4

Ref. no. 1895.1



ISO module trolley 7

Ref. no. 1899.1

Trolleys with metal drawers



Anaesthetic trolley 2

Ref. no. 1860.1



Treatment trolley 5

Ref. no. 1897.1



Treatment trolley 2

Ref. no. 1838.1

ISO module and function trolleys



ISO emergency trolley 2

Ref. no. 1858.1



ISO treatment trolley 3

Ref. no. 1864.1



ISO treatment trolley 8

Ref. no. 1761.1

Special models



Bronchoscopy trolley

yellow green RAL 6018 10785950



Treatment trolley

zinc yellow RAL 1018 10785700



Injection trolley

sky blue RAL 5015 10785800


Customised is our standard

For maximum mobility in the workplace

Since 1933, our company has been employing smart ideas and applying practical solutions to develop countless innovative product lines in the area of multipurpose equipment carts. These have been marketed worldwide with great success.

With our products and solutions, we guarantee maximum mobility, flexibility and reliability in the workplace. This means everyone can focus fully and completely on their own work.

We consistently implement the basic modular idea to develop practical solutions. The wide variety of standard components makes each individual cart an individually tailored customer solution – whether it's for a doctor's practice, small volume production or large quantities for equipment manufacturers. Even if you need to have it tailored just for you – HAEBERLE will make this possible.



HAEBERLE HAS IT: a modular basic system for your multifunctional carts. And if you need something extra? No problem, remember HAEBERLE MAKES IT: exactly the customised solution you require.

Variety of colours

Choose from 5 fashionable colours for decorative strips and handles. In your order please specify the colour code by adding it as the last digit to the reference number.



trolley/grey white
andles and decorative strips grey

Ref. no. xxxx.1



trolley/grey white
andles and decorative strips blue

Ref. no. xxxx.2



trolley/grey white
andles and decorative strips red violet

Ref. no. xxxx.3



trolley/grey white
andles and decorative strips yellow

Ref. no. xxxx.4



trolley/grey white
andles and decorative strips turquoise

Ref. no. xxxx.5

RAL colours

Or select your favourite colour for your keo® at an additional price.



System parts and accessories

HAEBERLE

mobile workstations



ISO module

Ref. no. 18897



Infusion bottle holder

Ref. no. 18020



Waste collector „tomo“

Ref. no. 18470



Anti-tamper seal for 3 grids

Ref. no. 17470



Catheter container

Ref. no. 18065



Stainless-steel holder for
dispenser boxes, for 1 box

Ref. no. 18707



ISO module
medication dispenser

Ref. no. 18400



Soap and disinfectant
dispenser, capacity 1litre

Ref. no. 17850



Organisation system with
2 inclined storage levels in
stainless steel

Ref. no. 17836



Disposal „Medi-Müll“ box
set with 20 boxes (1,5 ltrs.)
and mount

Ref. no. 18050



Waste collector „tango“
8 ltrs. with knee operated
opening system

Ref. no. 18358



Standard profile rail,
stainless steel for lateral
fastening to keo®

Ref. no. 17943

Additional system components and accessories can be found in the keo® price list.

Your medical distributor

HAEBERLE GmbH + Co. KG

Breitwiesenstr. 13, DE-70565 Stuttgart

Postfach 800524, DE-70505 Stuttgart

Telefon +49 (0) 711 78314-0

Telefax +49 (0) 711 78314-77

info@haeberle-med.de

www.haerberle-med.de