

HAEBERLE

mobile workstations

Type list valid from 1 July 2019








practical and ergonomical

furnishing accessories



PicBox[®] + Halux[®] + Mediskop[®]

Dispenser systems

Item	Description	Ref. no.
<h3>Dispenser systems PicBox® series</h3> <p>The tried and tested dispensers with the first-in / first-out system and dust protection according to DIN 58953. Suitable for syringes, canulas, dressing material or other small items. The boxes are filled from above, the material is removed from below. This guarantees that the oldest material is removed at first.</p> <p>The PicBox® series consists of grey-white (RAL 9002) plastic casing and dust-proof, transparent boxes.</p>		
Injection set "PicBox®"	 <p>Dispenser for syringes with 3 big boxes, dispenser for canulas with 5 small boxes.</p> <p>Dim.: 600 x 220 x 410 mm (w x d x h)</p> <p>Wall model</p>	11384
Injection set "PicBox® Plus"	 <p>Dispenser for syringes with 3 big boxes and 5 small boxes for canulas and with KOMBI set (with swab dispenser, swab moistener and additional container).</p> <p>Dim.: 725 x 235 x 410 mm (w x d x h)</p> <p>Wall model</p>	11380
Dispenser for syringes "PicBox®"	 <p>with 3 large big boxes.</p> <p>Dim.: 600 x 220 x 245 mm (w x d x h)</p> <p>Wall model</p>	11368
Dispenser for canulas "PicBox®"	 <p>with 5 small boxes.</p> <p>Dim.: 600 x 175 x 165 mm (w x d x h)</p> <p>Wall model</p>	11369
Injection set "PicBox® multi"	 <p>description see injection set PicBox®, but with white plastic veneered casing. Wall mounting rail for easy mounting included.</p> <p>Dim.: 630 x 235 x 435 mm (w x d x h)</p>	11481
Spare box big	 <p>suitable for syringe dispenser "PicBox®" with 3 big boxes, dust protection lid included.</p> <p>Dim.: 190 x 195 x 225 mm (w x d x h)</p>	11385
Spare box small	 <p>suitable for canula dispenser "PicBox®" with 5 small boxes, dust protection lid included.</p> <p>Dim.: 110 x 145 x 150 mm (w x d x h)</p>	11386
Spare dust protection lid	<p>for big spare box (ref. no. 11385) for small spare box (ref. no. 11386)</p>	14022 14023

Item	Description	Ref. no.
------	-------------	----------

Dispenser for syringes and canulas with simple boxes

Filling and removal from above. Made of impact-resistant polystyrene, recyclable and food-safe. Plastic casing, grey-white with transparent boxes. Temperature resistant from -20° C to +80° C.

Dispenser 2 for syringes



with 2 big boxes.

Dim.: 600 x 246 x 353 mm (w x d x h)

Wall model

05016

Spare box

05117

Dim.: 268 x 220 x 243 mm (w x d x h)

Dispenser 3 for syringes



with 3 big boxes.

Dim.: 600 x 163 x 242 mm (w x d x h)

Wall model

17208

Spare box

16391

Dim.: 175 x 140 x 160 mm (w x d x h)

Dispenser 4 for syringes



with 4 medium-sized boxes.

Dim.: 600 x 137 x 207 mm (w x d x h)

Wall model

17210

Spare box

16392

Dim.: 128 x 128 x 133 mm (w x d x h)

Dispenser 5 for canulas



with 5 small boxes.

Dim.: 600 x 110 x 164 mm (w x d x h)

Wall model


17209

Spare box

16393

Dim.: 100 x 97 x 108 mm (w x d x h)

Dispenser 6 for canulas



with 6 small boxes.

Dim.: 600 x 75 x 113 mm (w x d x h)

Wall model


05017

Spare box

05107

Dim.: 80 x 63 x 70 mm (w x d x h)

Dispenser 9 for canulas



with 9 small boxes.

Dim.: 600 x 51 x 77 mm (w x d x h)

Wall model








05018

Spare box

05108

Dim.: 50 x 42 x 46 mm (w x d x h)

Dispenser systems

Item	Description	Ref. no.
<h3>Dispenser systems COMPACT series</h3>	<p>Elegant dispensers for syringes, canulas or dressing material with ergonomic refill shafts and hygienic acrylic glass fronts. By uniform dimensions, the individual COMPACT elements can be optimally combined as contiguous wall installations. High quality cases in modern design. Ergonomically optimal arrangement of all components.</p> <p>The COMPACT series consists are all white casing, plastic veneered with dust-proof acrylic glass fronts. Dim.: 600 x 220 x 450 mm (w x d x h). Delivery includes wall mounting rails and mounting material.</p>	
Injection set "COMPACT"		<p>with practical refill shafts and inclined levels for syringes (3 compartments), canulas (5 compartments), such as 1 compartment for swap roll. Filling and removal from above - first-in/first-out.</p> <p style="text-align: right;">04877</p>
Injection set "COMPACT" 4 syringe compartments		<p>with practical refill shafts and inclined levels for syringes (4 compartments), canulas (4 compartments). Filling and removal from above - first-in/first-out.</p> <p style="text-align: right;">04879</p>
Dressing dispenser "COMPACT"		<p>with optimal arrangement for all common tube bandages, gauze bandages, plasters, etc. With 2 s/s holders for scissors, etc. One of them with magnetic support surface.</p> <p style="text-align: right;">04878</p>
Wall cupboard "COMPACT"		<p>with 1 intermediate compartment and 2 transparent acrylic glass doors.</p> <p style="text-align: right;">04876</p>
Organization cabinet "COMPACT"		<p>for forms with various dividers, without acrylic glass cover.</p> <p style="text-align: right;">14858</p>
Storage console "COMPACT"		<p>with storage compartment and telescope-guided drawer, incl. screws and dowels.</p> <p>Dim.: 600 x 450 x 280 mm (w x d x h)</p> <p style="text-align: right;">04872</p>
Wall desk "COMPACT"		<p>with inclined writing surface with insert for writing utensils. Removable pen holder. Storage compartment, 200 mm clear height. With wall support rail incl. screws and dowels.</p> <p>Dim.: 560 x 450 x 360 mm (w x d x h)</p> <p style="text-align: right;">15578</p>

Item	Description	Ref. no.
------	-------------	----------

Disposal boxes "Schluckfix" and "Medi-Müll"

For the safe, injury-free disposal of used canulas, ampoules, scalpels etc. The disposal boxes are made of plastic which is ecologically friendly and can be burnt without residues.

Container for canulas "Schluckfix"



Disposable container for approx. 500 used canulas. Made of black recycling plastic. Bulk discounts on request.

Dim.: 55 x 100 x 170 mm (w x d x h)

Set with 20 boxes
and table or wall mount

09506

set with 20 boxes
without mount

15502

set with 100 boxes
without mount

09490

Mount

separately for "Schluckfix"
for table or wall

04667

Disposal boxes "Medi-Müll" à 1.5 litres



made of yellow polypropylene with green lid. Offers 3 practical ways to dispose all types of canulas, vacutainers, scalpels, pipettes etc.

Dim.: Height 152 mm
Ø lid 143 mm
Ø base 108 mm

The boxes are delivered space-savily stacked inside each other, lid separately.

Set with 20 boxes
and wall mount

04780

set with 20 boxes
without mount

15501

set with 100 boxes
without mount

16653

Disposal boxes "Medi-Müll" à 2.5 litres



with an extremely big slot, which can be enlarged step by step by breaking out the panels. Ideally suited for scalpels, pipettes, broken glass and small dressing materials. Offers three practical ways for tearing-off all types of canulas, even for insulin canulas.




Dim.: Height 242 mm
Ø of lid 143 mm
Ø of base 108 mm

The boxes are delivered space-savily stacked inside each other, lid separately.



Original carton **with 45 boxes**, without mount







16654

Disposal boxes / Various dispenser items


Item	Description	Ref. no.
Holder	 <p>for "Medi-Müll" disposal boxes with 1,5 or 2,5 litres capacity, white-aluminium (RAL 9006), powder-coated.</p> <p>Wall model</p>	04671
Holder to standard profile rail	 <p>for "Medi-Müll" disposal box and for mounting to the right or left of all standard profile rails.</p> <p>Please order standard profile rail separately.</p>	16677
Disposal boxes "Medi-Müll" à 0.5 litre	 <p>Dim.: Height 110 mm Ø of lid 98 mm Ø of base 79 mm</p> <p>Boxes are delivered space-savingly stacked inside of each other, lid separately.</p> <p>Original carton with 170 boxes</p>	16655
Adhesive adapter	<p>suitable for "Medi-Müll" boxes of 0.5/ 1.5/2.5 litres capacity; for fixing the boxes on all surfaces, such as tables or trays.</p> <p>Minimum order quantity: 10 pieces</p>	17212

Various dispenser items



Soap and disinfectant dispenser	 <p>made of naturally anodised aluminium, for standard bottles with 1 litre capacity.</p> <p>Dim.: 95 x 185 x 290 mm (w x d x h)</p>	17851
Drip tray	<p>separately, removable; to clip on soap and disinfectant dispenser (ref. no. 17851).</p> <p>Dim.: 150 x 260 x 410 mm (w x d x h)</p>	16256
Sensor dispenser	 <p>with hygienic sensor technology for the contact-free dispensing of liquid soap and disinfectants. Made of naturally anodized aluminium, for standard bottles with 1 litre capacity.</p> <p>For wall mounting with fixing material</p> <p>Dim.: 95 x 210 x 290 mm (w x d x h)</p>	11140
Drip tray	<p>separately, removable; to clip on sensor dispenser (ref. no. 11140).</p> <p>Dim.: 150 x 260 x 410 mm (w x d x h)</p>	16256

Item	Description	Ref. no.
KOMBI set	 <p>for blood taking and injection; consisting of swab dispenser "TAKEFIX", swab moistener "TIPPFIX" and additional container. All parts are made of unbreakable, grey-white plastic material.</p> <p>With mounting material for fastening on syringes dispenser "PicBox®" (ref. no. 11368).</p> <p>Dim.: 120 x 230 x 265 mm (w x d x h)</p>	04767
		16580
Swab dispenser "TAKEFIX"	 <p>with casing made of unbreakable, grey-white (RAL 9002) plastic.</p> <p>Dim.: 120 x 230 x 205 mm (w x d x h)</p>	04803
Swab moistener "TIPPFIX"	 <p>with casing made of unbreakable, grey-white plastic (RAL 9006). Filling capacity 250 ccm.</p> <p>Dim.: 60 x 150 x 60 mm (w x d x h)</p> <p>Bulk discounts on request.</p>	09493
Ampoule opener "CUPFIX"	 <p>Receptacle for used ampoules, removable. Diamond file is springloaded in a stainless steel plate. Plastic housing, grey-white,</p> <p>Dim.: 70 x 150 x 200 mm (w x d x h)</p> <p>Wall model, incl. assembly set</p> <p>Spare diamond file for "CUPFIX"</p>	09507
		04663
Holder for dispenser boxes	 <p>made of brushed stainless steel. Suitable for 1 to 4 common dispenser boxes of disposable gloves or tissues. For wall mounting (without mounting material).</p> <p>Holder 1-fold (for 1 box) Inner dim.: 142 x 94 x 250 mm (w x d x h)</p> <p>Holder 3-fold (for 3 boxes) Inner dim.: 253 x 94 x 394 mm (w x d x h)</p> <p>Holder 4-fold (for 4 boxes) Inner dim.: 253 x 94 x 525 mm (w x d x h)</p>	11222
		11223
		05074
Universal holder "Snap"	 <p>made of white plastic, for dispenser boxes of all kinds. Ideal for boxes with disposable gloves, tissues or wipes.</p> <p>The spring-loaded clamp is adjustable to fit many box sizes up to 1 kg. Fixing with double adhesive stripes.</p> <p>Dim.: 100 x 20 x 165 mm (w x d x h)</p> <p>With individual imprint ideal as a promotional gift. Bulk discounts on request.</p>	15572

ENT products

Item	Description	Ref. no.
Tray for medications "Langenau"	 <p>including 1 prescription pad. Base made of white and cover made of transparent plastic, 18 mm high.</p> <p>with 16 compartments dim.: 310 x 310 mm (w x d)</p>	04705
	<p>with 20 compartments dim.: 310 x 390 mm (w x d)</p>	04706
	<p>with 25 compartments dim.: 390 x 390 mm (w x d)</p>	04707
Spare cover	<p>made of transparent plastic for "Langenau" tray for medications,</p> <p>16 compartments</p>	05247
	<p>20 compartments</p>	05248
	<p>25 compartments</p>	05249
Spare prescription pad	<p>for "Langenau" tray for medications. 50 sheets with 300 tear off slips. Per pad</p> <p>Minimum order quantity: 20 pads.</p>	04708

ENT products

Mirror heater "Mirrotherm® Plus"	 <p>with extremely long service life of the heating elements.</p> <p>An additional blower ensures fast heating up and automatic cooling down. Easy operation. Luminous indication (LEDs) GS-marked. Plastic case, umbergrey.</p> <p>Dim.: 213 x 120 x 95 mm (w x d x h)</p>	08131
Dispenser box for "KREPPFIX" tongue swabs	 <p>The dispenser box consists of a white bottom in sturdy plastic with pivot for the fixation of the "KREPPFIX" tongue swab roll and a transparent top with opening slit.</p> <p>The box may be fixed with 3 suction cups even on smooth surfaces. Height 10 cm, Ø 12.5 cm. Delivery without "KREPPFIX" tongue swabs.</p>	04711
Tongue swabs "KREPPFIX"	<p>Crepe paper strips, natural white, perforated, approx. 60 mm wide, Ø 90 mm. 25 m long, approx. 150 sheets. Length per single sheet approx. 160 mm. Central pharmaceutical no. 10299603. Minimum order quantity: 3 packages.</p> <p>Package with 6 rolls</p> <p>Original carton with 18 packages of 6 rolls</p>	04709 16395
Absorbent cotton	<p>with removal slot. Length 4 m, Ø 95 mm, approx. 60 mm wide, rolled. Minimum order quantity: 6 rolls.</p> <p>Per roll</p>	04972

Item	Description	Ref. no.
------	-------------	----------

Medical lamps

Whether for doctor's practices, clinics, nursing homes or cosmetics studios - The high quality lamps are suitable for many application areas. All lamps comply with the latest version of medical norms and are CE tested. EN 60601, EN 60598, EN 12464. The reference numbers do not include mounts, unless mentioned otherwise. Suitable mounts can be found on page 10.

Halux 35 F ENT back-mirror lamp



The practical shapely **halogen back-mirror lamp with frosted glass filter** and spring-balanced articulated arm.

Lamp:	1 x halogen bulb
Operating unit:	electronic transformer
Supply voltage:	230 V/50 Hz
Power supply cable:	1.4 m Euro plug
Luminaire housing:	polyamid (PA)
Weight (net):	1,3 kg
Power consumption:	40 W
Switch:	two-step rocker switch
Norms:	EN 60601-1, EN 60601-2-41
Joint balance:	spring
Special feature:	frosted glass filter
Protection class:	II

16425

LED examination lamps

The use of LED technology in medical lamps provides high illumination intensity at low heat release, long service life, high efficiency and excellent quality of light.

The LEDs have a minimum service life of 30.000 operating hours.

3 colour shades:	3400 K (warm white)	- for different types of skin examination
	3800 K (neutral white)	- for general examinations
	4300 K (cold light)	- ideal for the recognition of different vessels, for post operative treatment and small interventions

The colour temperature - measured in Kelvin (K) gives information about the light colour.

Halux® LED F 7-P





The robust **LED examination lamp** meets all important criteria for precise examinations. The spring balanced swivel arm enables quick positioning. The lamp head has a double-wall casing for optimized cooling.

Lamp:	1 x LED module
Light output:	15.000 Lux/0.5 m
Illumination field:	Ø 30 cm / 0.5 m
Colour temperatur:	4000 K
Operating unit:	electronic plug-in transformer
Supply voltage:	100 - 240 V; 50/60 Hz
Power supply cable:	2.0 m Euro plug
Luminaire housing:	polyamide (PA)
Weight (net):	approx. 1,7 kg
Power consumption:	19 - 27 VA
Switch:	rocker switch (I/O)
Joint balance:	friction
Norms:	EN 60601-1, EN 60601-2-41
Special feature:	emission angle 10°
Protection class:	II

05522

Med. lamps

Item	Description	Ref. no.
Visiano LED S	 <p>The LED examination lamp has a slimline design with ergonomic handle and a bright, homogeneous light thanks to unique diamond optics. This lamp convinces with its excellent illumination intensity and a large light field. The lamp head can be rotated around 280° and has a closed casing with a smooth surface for easy cleaning.</p> <p>Lamp: 1 x LED module Light output: 50.000 Lux/0.5 m Illumination field: Ø 18 cm/0.5 m Colour temperatur: 4400 K Operating unit: electronic converter Supply voltage: 100 - 240 V; 50/60 Hz Power supply cable: 1.8 m Euro plug Luminaire housing: PMMA Weight (net): approx. 1,1 kg Power consumption: 13 W Switch: rocker switch (I/O) at the arm Norms: EN 60601-1, EN 60601-2-41 Joint balance: friction Special feature: with handle Protection class: II</p>	05521

Culta LED S	 <p>The glare-free light and extra large light field provide maximum reading comfort with a pleasant light colour. The LED reading lamp has a satin-transparent lampshade and harmonizes with any interior design. The lamp head can be rotated around 240°. With the gooseneck the positioning is easy and flexible.</p> <p>Lamp: LED module Illumination field: Ø 60 cm/0.5 m Colour temperatur: 3000 K Operating unit: electronic plug-in transformer Supply voltage: 100 - 240 V; 50/60 Hz Power supply cable: 2.0 m Luminaire housing: polycarbonat (PC), satine Weight (net): 1,0 kg Power consumption: 5,2 W Switch: cord switch Norms: EN 60598-1, EN 60598-2-25 Special feature: Gooseneck white-aluminium with satin-transparent lampshade Protection class: II</p>	05523
--------------------	--	--------------

Mounts for table / wall / standard profile rail

Table holder	 <p>with removal guard made of polyamid GF, white (RAL 9010). Suitable for lamps with stud Ø 16 mm. Weight: 0,09 kg Dim.: 60 x 60 x 60 mm (w x d x h)</p>	09947
Wall holder	 <p>with removal guard and cable hole, made of polyamid GF, pure white (RAL 9010). Suitable for lamps with stud Ø 16 mm, weight: 0,14 kg.</p>	15899
Rail bracket	 <p>with removal guard for standard profile rail, made of eloxal aluminium, suitable for lamps with stud Ø 16 mm. Weight: 0,16 kg.</p>	09949

Mediskop®- X-ray film viewers

An eye-catcher for your medical office is offered by the **Mediskop® device types 560** and **860** with on/off-switch or the model **560-HF** with luminosity regulation and high frequency technology, which offers up to 50% longer lifetime of the fluorescent lamps during operation.



Silver-grey quality aluminium housing with grey salient angles • Front panel made of acrylic glass • All-around film holding rail for maximum viewing comfort • Supply voltage 230 V / 50 Hz • Lowering of the luminance down to 500 cd/m² • Screwless front for installation and service easily removable.

Mediskop® with on/off switch

Type	Viewing surface w x h (mm)	External dim. w x d x h (mm)	Number of tubes (Watt)	Wattage (Watt)	Colour temp. of transmission light (Kelvin)	Colour temperature of tubes (Kelvin)	Illuminance max. (Lux)	luminance (cd/qm approx.)	Weight (kg approx.)	Colour of angles	Ref. no.
560	560 x 435	635 x 100 x 510	2/18	45	6000	6500	6200	1790	9,4	grey	15221
860	865 x 435	940 x 100 x 510	2/30	80	6000	6500	9300	2630	12,8	grey	15223

Mediskop® with stepless luminosity regulation and high frequency technology

Type	Viewing surface w x h (mm)	External dim. w x d x h (mm)	Number of tubes (Watt)	Wattage (Watt)	Colour temp. of transmission light (Kelvin)	Colour temperature of tubes (Kelvin)	Illuminance max. (Lux)	luminance (cd/qm approx.)	Weight/ (kg approx.)	Colour of angles	Ref. no.
560 HF	560 x 435	635 x 100 x 510	2/18	35	6000	6500	7350	2100	9,3	grey	17300

All prices are understood exclusive of VAT. As an innovative company we reserve the right to amend prices and designs. As for the rest our General Terms of Business in the valid version apply.


WEEE-Reg. no. DE 20235235



HAEBERLE

mobile workstations

Our versatile product range



Your medical distributor

HAEBERLE GmbH + Co. KG

Breitwiesenstr. 13, DE-70565 Stuttgart

P.O.Box 800524, DE-70505 Stuttgart

Phone +49 (0) 711/78314-0

Fax +49 (0) 711/78314-77

info@haeberle-med.de

www.haerberle-med.de

Version 07/19