

HAEBERLE

mobile workstations

practical and ergonomical

furnishing accessories



PicBox[®] + Halux[®] + Mediskop[®]

Furnishing accessories, practical and ergonomical

PicBox® series

The tried and tested dispensers with the first-in / first-out system and dust protection according to DIN 58953. Suitable for syringes, canulas, dressing material or other small items.

COMPACT series

Elegant dispensers for syringes, canulas or dressing material with ergonomic inclined planes and hygienic acrylic glass fronts. The uniform dimensions of the COMPACT elements means they can be combined optimally to create large wall collections.

Disposal boxes and various dispensers

Disposal articles „Medi-Müll“ and „Schluckfix“ as well as practical auxiliaries for taking blood samples, disinfection and trays for medications.

ENT products

Include the popular mirror heater „Mirrotherm® Plus“, tongue swabs with matching dispenser box and the halogen back-mirror lamp „Halux 35 F“.

Medical lamps

High-quality, modern LED lamps for examination and treatment. All lamps comply with the hygienic requirements and the relevant medical standards.

X-ray film viewers

„Mediskop“ devices with elegant quality aluminium housing in 2 sizes and grey salient angles. Available with on/off-switch or with luminosity regulation and high frequency technology (HF).




The features

- PicBox® series with grey-white plastic casing and transparent dispensing boxes with practical removal opening.
- First-in / first-out system.
The boxes are filled from above and the material is taken from below.
- COMPACT dispensers elegant and white with hygienic acrylic glass fronts.
- Medi-Müll disposal boxes for the injury-free disposal of used canulas, scalpels, etc.
- Lamps with optimum illumination and high-contrast light due to LED technology.
- Individual illumination level, long service life for the LED's, low running costs.
- The X-ray viewers offer an optimal, even illumination and a large brightness control range as well as fatigue-free working with high frequency technology.

Dispenser systems PicBox® series

Injection set PicBox® / PicBox® Plus / PicBox® multi



Injection set PicBox®

Ref. no. 11384



Injection set PicBox® Plus

Ref. no. 11380



Injection set PicBox® multi

Ref. no. 11481

Syringe dispenser PicBox® / Canula dispenser PicBox® / Spare boxes



Dispenser for syringes PicBox®

Ref. no. 11368



Dispenser for canulas PicBox®

Ref. no. 11369



Spare boxes PicBox®

Ref. no. big / small 11385 / 11386

Dispenser for syringes / Dispenser for canulas with simple boxes



Dispenser for syringes 2

Ref. no. 05016



Dispenser for syringes 4

Ref. no. 17210



Dispenser for canulas 6

Ref. no. big / small 05017

and COMPACT series

COMPACT dressing dispenser / injection sets



COMPACT dressing dispenser

Ref. no. 04878



COMPACT injection set


Ref. no. 04877



COMPACT injection set
with 4 syringe compartments

Ref. no. 04879

COMPACT organizer cupboard / wall cupboard / wall console



COMPACT organizer cupboard

Ref. no. 14858



COMPACT wall cupboard

Ref. no. 04876



COMPACT wall console

Ref. no. 04872

LED - examination lamps / reading lamp



Visiano LED S

Ref. no. 05521



Halux® LED F 7-P

Ref. no. 05522



Culta LED S

Ref. no. 05523

Mediskop®

X-ray film viewers

HAEBERLE

An eye-catcher for your medical office is offered by the **Mediskop®-device types 560 und 860** with on/off-switch or the model **560-HF** with luminosity regulation and high frequency technology.

- The designkick for doctor's office or clinic
- Optimal, uniform illumination
- Available with or without luminosity regulation
- Even illumination and great luminosity range
- Strain-free working due to high frequency technology (HF)
- HF devices offer higher operating efficiency. Up to 50% longer lifetime of the fluorescent lamps during operation.



Features:

- Silver-grey quality aluminium housing with grey salient angles
- Front panel made of acrylic glass
- All-around film holding rail for maximum viewing comfort
- Supply voltage 230 V / 50 Hz
- Lowering of the luminance down to 500 cd/m²
- Screwless front for assembly and service easily removable

ENT products / Dispensers

Disposal boxes

HAEBERLE



Halux 35 F ENT back-mirror lamp

Ref. no. 16425



Mirror heater „Mirrotherm® Plus“

Ref. no. 08131



Dispenser box for tongue swabs

Ref. no. 04711



Swab moistener „TIPPFIX“

Ref. no. 09493



Swab dispenser „TAKEFIX“

Ref. no. 04803




KOMBI-set

Ref. no. 04767



Ampoule opener „CUPFIX“

Ref. no. 09507



Trays for medications „Langenau“

Ref. no. 04705
with 16 compartments

Ref. no. 04706
with 20 compartments

Ref. no. 04707
with 25 compartments



Holder for dispenser boxes

Ref. no. 11222
for 1 box

Ref. no. 11223
for 3 boxes

Ref. no. 05074
for 4 boxes



Soap and disinfectant dispenser

Ref. no. 11140
Sensor dispenser

Ref. no. 16256
Drip tray



Disposal boxes „Medi-Müll“

Ref. no. 16655
0,5 ltrs. / 170 boxes

Ref. no. 16653
1,5 ltrs. / 100 boxes

Ref. no. 16654
2,5 ltrs. / 45 boxes



Container for canulas „Schluckfix“

Ref. no. 09506
20 boxes with
wall- / table mounting

Ref. no. 15502
20 boxes
without holder

Ref. no. 09490
20 boxes
without holder

Further items can be found in the price list furnishing accessories.

Your medical distributor

HAEBERLE GmbH + Co. KG

Breitwiesenstr. 13, DE-70565 Stuttgart

P.O.Box 800524, DE-70505 Stuttgart

Phone +49 (0) 711/78314-0

Fax +49 (0) 711/78314-77

info@haeberle-med.de

www.haerberle-med.de