

HAEBERLE

mobile workstations

doppio
by HAEBERLE

puristical and
powerful

equipment cart system

type list

valid from 1. October 2019

doppio

- the puristical and powerful cart system

The equipment cart system **doppio** offers you an individual design option for the doctor's practice or clinic.

The modern, multifunctional cart system convinces with high-quality, coated shelves made of wood/plastic with rounded edges. The shelves are available with widths of 45 cm and 60 cm and offer enough space for any requirements.

The newly designed aluminium carrier profiles with integrated cable ducts and coloured decorative strips in grey or blue, make the carts the ideal mobile complement to your practice equipment.

The highly sophisticated design guarantees aesthetic integration of the cart systems into modern practice furniture. The large selection of system parts and accessories, allow plenty of possibilities for individual design.

**HAEBERLE
HAS IT**

**HAEBERLE
MAKES IT**

Features: Powerful chassis, 4 elegant aluminium carrier profiles with cable ducts and for attaching accessories. High-quality, beautifully shaped, white shelves and drawers with rounded soft edges.

Function: As basic cart, equipment cart, ergometry cart, tablet/iPad cart, injection cart ... and more.

Application: In clinics, rehabilitation centres and doctors' practices.



POAG equipotential bonding system conforming to DIN 42801

In rooms used for medical purposes, equipment and other conductive parts may have different potentials, which poses a risk to persons. Equipotential bonding offers safe protection.

Optionally, POAG pins will be mounted ex works on equipment carts and electrically tested during the final inspection. If your rooms do not already have one, please order the correct equipotential bonding cable as well. Our team would be happy to advise you.

The POAG pins will be attached at rear of a horizontal arm and between the carrier profiles of the equipment carts. Retrofitting by the customer is possible, but must be carried out by qualified personnel.

Components

- Sturdy chassis with large quality light-running castors.
- Electrically conductive castors optionally for an extra charge.
- Puristical aluminium carrier profiles with profile cover caps and integrated cable ducts for convenient cable management.
- Decorative strips in the profiles in grey or blue available.
- Loadable, white-coated shelves, made of wood/plastic with rounded edges.
- Beautifully shaped drawers with rounded edges, grey drawer handle and easy-sliding ball telescopic guide with snap-in mechanism.
- Handy, lateral push handles in grey.

Technical details / dimensions

Cart outer dimensions:	With shelf width 450 mm and castors Ø 100 mm 565 x 571 x 793 mm (w x d x h) 565 x 571 x 1177 mm (w x d x h) 565 x 571 x 1474 mm (w x d x h) with shelf width 600 mm and castors Ø 100 mm 715 x 571 x 793 mm (w x d x h) 715 x 571 x 1177 mm (w x d x h) 715 x 571 x 1474 mm (w x d x h)
Maximum load:	up to 120 kg
Carrier profiles height:	566 / 950 / 1206 mm as standard
Shelf outer dimensions:	450 x 450 x 22 mm (w x d x h) 600 x 450 x 22 mm (w x d x h)
Shelf load limit:	up to 30 kg per shelf
Drawer outer dimensions:	450 x 450 x 100 mm (w x d x h) 600 x 450 x 100 mm (w x d x h)
Drawer surface load:	up to 10 kg
Quality light-running castors:	4 castors Ø 100 mm or Ø 125 mm with wall-guards as standard. Castors with Ø 150 mm and electrically conductive castors are available for an extra charge.

Materials / colours

Carrier profiles:	Anodized aluminium
Cover cap:	Polyamide, grey (6502B acc. NCS)
Chassis:	Aluminium cast, fine textured, powder-coated, grey (RAL 7012)
Decorative strips:	Made of PVC, grey (6502B acc. NCS) or blue (3502B acc. NCS)
Shelves:	Wood/plastic, coated, white (W980-ST2 EGGER)
Push handle:	Polyamide, grey (6502B acc. NCS)
Drawer handle:	Polyamide, grey (6502B acc. NCS)
Drawer:	Wood/plastic, coated, white (W980-ST2 EGGER)

For all cart models listed in the price list, the illustrated push handles are exclusive. All listed carts are delivered assembled at no extra charge and can be retrofitted at any time.

Basic cart / Equipment carts

doppio - basic cart



with 2 shelves, 1 drawer with drawer handle, 4 castors \varnothing 100 mm.
Up to 4 drawers can be retrofitted.

Cart 45

Dim.: 565 x 571 x 793 mm (w x d x h)

Decorative strips	Ref. no.
grey	26101
blue	26102

Cart 60

Dim.: 715 x 571 x 793 mm (w x d x h)

Decorative strips	Ref. no.
grey	26201
blue	26202

Single drawers see on page 9.

doppio - equipment cart 1



with 3 shelves, 1 horizontal arm, 4 castors \varnothing 100 mm with 2 locking brakes.

Cart 45

Dim.: 565 x 571 x 1177 mm (w x d x h)

Decorative strips	Ref. no.
grey	26341
blue	26342

Cart 60

Dim.: 715 x 571 x 1177 mm (w x d x h)

Decorative strips	Ref. no.
grey	26351
blue	26352

doppio - equipment cart 2



with 4 shelves, 1 horizontal arm, 1 drawer with drawer handle, 4 castors \varnothing 100 mm with 2 locking brakes.

Cart 45

Dim.: 565 x 571 x 1177 mm (w x d x h)

Decorative strips	Ref. no.
grey	26211
blue	26212

Cart 60

Dim.: 715 x 571 x 1177 mm (w x d x h)

Decorative strips	Ref. no.
grey	26111
blue	26112

Endoscopy cart / Ergometry cart / Tablet/iPad cart

doppio - endoscopy cart



with 5 shelves, 1 horizontal arm, 1 drawer with drawer handle, 4 electrically conductive castors (cart 45 / Ø 125 mm, cart 60 / Ø 150 mm) with 2 locking brakes.

Cart 45

Dim.: 565 x 571 x 1474 mm (w x d x h)

Decorative strips	Ref. no.
grey	26221
blue	26222

Cart 60

Dim.: 715 x 571 x 1504 mm (w x d x h)

Decorative strips	Ref. no.
grey	26241
blue	26242

doppio - ergometry cart



with 2 shelves (450 mm deep), 1 shelf at top (350 mm deep), 1 drawer with drawer handle, 1 cable and hose holder (ref. no. 13926) incl. attachment device, 4 castors Ø 100 mm with 2 locking brakes.

Cart 45

Dim.: 565 x 571 x 1177 mm (w x d x h)

Decorative strips	Ref. no.
grey	26231
blue	26232

Cart 60

Dim.: 715 x 571 x 1177 mm (w x d x h)

Decorative strips	Ref. no.
grey	26301
blue	26302

doppio - tablet/iPad cart



with 3 shelves, 1 horizontal arm, 1 drawer with drawer handle, 1 mount for tablet/iPad (ref. no. 20854), 4 castors Ø 100 mm with 2 locking brakes.

Cart 45

Dim.: 565 x 571 x 1177 mm (w x d x h)

Decorative strips	Ref. no.
grey	26321
blue	26322

Cart 60

Dim.: 715 x 571 x 1177 mm (w x d x h)

Decorative strips	Ref. no.
grey	26331
blue	26332

Anaesthesia cart / Injection carts

doppio - anaesthesia cart



with 2 shelves, 5 drawers with drawer handles, 4 electrically conductive castors Ø 125 mm with 2 locking brakes.

Cart 45

Dim.: 565 x 571 x 834 mm (w x d x h)

Decorative strips	Ref. no.
grey	26311
blue	26312

Cart 60

Dim.: 715 x 571 x 834 mm (w x d x h)

Decorative strips	Ref. no.
grey	26271
blue	26272

doppio - injection cart COMPACT



with 2 shelves, 1 drawer with drawer handle, 1 drawer inset with variable partitions, 1 injection set COMPACT, disposal boxes "Medi-Müll" 1.5 ltrs. (set with 20 boxes) with mount, 1 waste container with s/s bowl and 1 arm support, each with attachment device, 4 castors Ø 100 mm with 2 locking brakes.

Cart 60

Dim.: 715 x 571 x 1433 mm (w x d x h)

Decorative strips	Ref. no.
grey	26261
blue	26262

doppio - injection cart PicBox-Plus



with 2 shelves, 1 drawer with drawer handle, 1 drawer inset with variable partitions, 1 injection set PicBox with KOMBI set, 1 canula container "Schluckfix", (set with 20 boxes), 1 waste container with s/s bowl and 1 arm support, each with attachment device, 4 castors Ø 100 mm with 2 locking brakes.

Cart 60

Dim.: 715 x 571 x 1433 mm (w x d x h)

Decorative strips	Ref. no.
grey	26251
blue	26252

System parts

Chassis / Carrier profiles

Chassis doppio

with basic shelf, 4 quality light-running castors with wall-guards.

for **cart 45**

Dim.: 565 x 571 x 223 mm (w x d x h)

Chassis	Ref. no.
Castors Ø 100 mm	26036

Dim.: 565 x 571 x 223 mm (w x d x h)

Chassis	Ref. no.
Castors Ø 100 mm 2 lockable	26039

Dim.: 565 x 571 x 265 mm (w x d x h)

Chassis	Ref. no.
Castors Ø 125 mm 2 lockable	26045

Dim.: 565 x 571 x 265 mm (w x d x h)

Chassis	Ref. no.
Castors Ø 125 mm electr. conductive, 2 lockable	26083

for **cart 60**

Dim.: 715 x 571 x 223 mm (w x d x h)

Chassis	Ref. no.
Castors Ø 100 mm	26035

Dim.: 715 x 571 x 223 mm (w x d x h)

Chassis	Ref. no.
Castors Ø 100 mm 2 lockable	26023

Dim.: 715 x 571 x 265 mm (w x d x h)

Chassis	Ref. no.
Castors Ø 125 mm 2 lockable	26043

Dim.: 715 x 571 x 265 mm (w x d x h)

Chassis	Ref. no.
Castors Ø 125 mm electr. conductive, 2 lockable	26084


Dim.: 715 x 571 x 294 mm (w x d x h)

Chassis	Ref. no.
Castors Ø 150 mm electr. conductive, 2 lockable	26044



Carrier profiles in standard heights

2 pairs of carrier profiles, made of aluminium, with integrated cable ducts, grey cover caps and decorative strips in grey or blue. For the guarantee of an optimal stability, each cart with carrier profiles up from 950 mm and without a drawer has to be equipped with one horizontal arm (ref. no. 26067) or (ref. no. 26068). See on page 9.



Profile height: 566 mm

Decorative strips	Ref. no.
grey	26021
blue	26022

Profile height: 1206 mm

Decorative strips	Ref. no.
grey	26081
blue	26082

Profile height: 950 mm


Decorative strips	Ref. no.
grey	26061
blue	26062

Profile height: 566 / 1206 mm

Decorative strips	Ref. no.
grey	26121
blue	26122

Carrier profiles in individual heights

4 carrier profiles made of aluminium, with integrated cable ducts, grey cover caps and decorative strips in grey or blue.



height \ colours	Decorative strips: grey	Decorative strips: blue
0 - 566 mm	26131	26132
566 - 1206 mm	26141	26142
1206 - 1846 mm	26151	26152

Castors / Push handle / Shelves / Anti-slide mat

Castors Ø 100 mm



Extra charge for 2 lockable plastic double swivel castors.

Ref. no.

13906

Castors Ø 100 mm electrically conductive

Extra charge for 4 plastic double swivel castors, 2 of them are lockable.

Ref. no.

13907

Castors Ø 125 mm



Extra charge for 4 plastic castors, 2 of them are lockable.

Ref. no.

26074


Castors Ø 125 mm electrically conductive

Extra charge for 4 plastic castors, 2 of them are lockable.

Ref. no.

26075

Castors Ø 150 mm



Extra charge for 4 metal castors, 2 of which lockable.

Ref. no.

26076

Castors Ø 150 mm electrically conductive

Extra charge for 4 metal castors, 2 of which lockable.

Ref. no.

26077

Push handle



grey, made of polyamide, 1 piece.

Ref. no.

26028

Shelf



350 mm deep, white, made of wood/plastic, with rounded edges.

for **cart 45**

Dim.:
450 x 350 x 22 mm (w x d x h)

Ref. no.

26032

for **cart 60**

Dim.:
600 x 350 x 22 mm (w x d x h)

Ref. no.

26025

Shelf



450 mm deep, white, made of wood/plastic, with rounded edges.

for **cart 45**

Dim.:
450 x 450 x 22 mm (w x d x h)

Ref. no.

26030

for **cart 60**

Dim.:
600 x 450 x 22 mm (w x d x h)

Ref. no.

26024

Anti-slide mat for shelves

made of ethylene vinyl acetate copolymer (EVA), transparent, noise-damming.

for **cart 45**

Dim.: 430 x 430 mm (w x d)

Ref. no.

26087

for **cart 60**

Dim.: 580 x 430 mm (w x d)

Ref. no.

26089

Dim.: 430 x 330 mm (w x d)

Ref. no.

26088

Dim.: 580 x 330 mm (w x d)

Ref. no.

26090

Drawer / Drawer inset / Traverse / Horizontal arms

Drawer

made of wood/plastic, white, with grey drawer handle, made of polyamide. Drawer surface load up to 10 kg.



for **cart 45**

Inner dim.:
376 x 400 x 80 mm (w x d x h)

Ref. no.

26034

for **cart 60**

Inner dim.:
526 x 400 x 80 mm (w x d x h)

Ref. no.

26026

Drawer inset

with variable partitions and flexible dividers. Made of transparent synthetic material, removable and suitable for cleaning in the dishwasher.



for **cart 45**

with 1 lengthways divider
and 4 crossways dividers
Dim.: 376 x 400 x 55 mm (w x d x h)

Ref. no.

26085

for **cart 60**

with 2 lengthways dividers
and 6 crossways dividers
Dim.: 526 x 400 x 55 mm (w x d x h)

Ref. no.

26086

Traverse wood/plastic

as additional rear upstand of the plate or tilt protection, white coated (W980-ST2 EGGER). **Not suitable for stabilisation.**



for **cart 45**

Dim.: 450 x 22 x 110 mm (w x d x h)

Ref. no.

26060

for **cart 60**

Dim.: 600 x 22 x 110 mm (w x d x h)

Ref. no.

26063

Horizontal arm metal

for additional stabilisation of the carts or for fastening a multiple socket outlet. Made of metal, powder-coated, grey-white (RAL 9002).



for **cart 45**

Dim.: 450 x 35 x 100 mm (w x d x h)

Ref. no.

26067

for **cart 60**

Dim.: 600 x 35 x 100 mm (w x d x h)

Ref. no.

26068

Horizontal arm metal

with 2 POAG pins (ref. no 26073). Up to 7 POAG pins are possible. Made of metal, powder-coated, grey-white (RAL 9002).



for **cart 45**

Dim.: 450 x 52 x 100 mm (w x d x h)

Ref. no.

26071

for **cart 60**

Dim.: 600 x 52 x 100 mm (w x d x h)

Ref. no.

26072

Horizontal arm metal

with 2 POAG pins (ref. no 26073) and a multiple socket outlet (ref. no 21868). Up to 7 POAG pins are possible. Made of metal, powder-coated, grey-white (RAL 9002).



for **cart 45**

Dim.: 450 x 79 x 100 mm (w x d x h)

Ref. no.

26078

for **cart 60**

Dim.: 600 x 79 x 100 mm (w x d x h)

Ref. no.

26079

"alegro" monitor carrier system - different variants

For stepless and precise fixation without any physical exertion. Horizontally and vertically rotatable, tilting and swivelling function. Monitor adaption VESA 75/100.

"alegro" monitor holder



on tightly adjusted rod in the carrier profile.
Height approx. 400 mm above the profile.

Load capacity up to 12 kg

Ref. no.

26125

Load capacity up to 18 kg

Ref. no.

26126

"alegro" holder for double monitoring



on aluminium-traverse (ref. no. 22010) between two carrier profiles, for optimal synchronous positioning of monitors. Cable duct inside the profiles.

for **cart 60**

Load capacity up to 2 x 12 kg

Ref. no.

25003

Load capacity up to 2 x 18 kg

Ref. no.

25029

"alegro" monitor holder



on aluminium-traverse (Ref. no. 22010) between two carrier profiles.

for **cart 45**

Load capacity up to 12 kg

Ref. no.

25064

Load capacity up to 18 kg

Ref. no.

25055

for **cart 60**

Load capacity up to 12 kg

Ref. no.

25062

Load capacity up to 18 kg

Ref. no.

25063

"alegro" monitor holder



on metal horizontal arm (ref. no. 26067 / 26068) on height-adjustable rod, between two carrier profiles.

for **cart 45**

Load capacity up to 12 kg

Ref. no.

26096

for **cart 60**

Load capacity up to 12 kg

Ref. no.

26097

"alegro" monitor holder on swivel arm



with 1 hinge joint and integrated cable duct. Reach 300 mm.

Load capacity up to 12 kg

Ref. no.

25077

Load capacity up to 18 kg

Ref. no.

25078

"alegro" monitor holder on swivel arm



with 2 hinge joints and integrated cable duct. Reach 300 mm.

Load capacity up to 12 kg

Ref. no.

25087

Load capacity up to 18 kg

Ref. no.

25088

"alegro" monitor holder on articulated arm



with adapter angle for fixing to the carrier profile and with integrated cable duct. Reach 300 mm.

Load capacity up to 12 kg

Ref. no.

25032

Load capacity up to 18 kg

Ref. no.

25035

"alegro" monitor holder on double articulated arm



with adapter angle for fixing to the carrier profile and with integrated cable duct. Reach 2 x 150 mm.

Load capacity up to 12 kg

Ref. no.

25036

Load capacity up to 18 kg

Ref. no.

25076

Accessories

TFT mounts / Swivelling arm with TFT mount / Waste container / Mount for tablet/iPad / Waste collectors

TFT mount



for flat screens (VESA 75/100), laterally attached to the profile. Rotatable and tiltable on height-adjustable rod. Steel, powder-coated, white-aluminium (RAL 9006). Load capacity up to 12 kg.

Ref. no.

26095

Swivelling arm with TFT mount



for flat screens (VESA 75/100), laterally attached to the profile. Rotatable and tiltable. Arm can be swivelled (range 300 mm), white-aluminium (RAL 9006). Load capacity up to 12 kg.

Ref. no.

26093

TFT mount



on height-adjustable rod for mounting on horizontal arm. Please order the horizontal arm (ref. no 26067 / 26068) separately.

Ref. no.

26123

Waste container



with s/s bowl, Ø 320 mm and attachment device for the profile.

Ref. no.

13911

Mount for tablet/iPad



for fixing to horizontal arm, only mountable ex works. Anodised aluminium, adjustable from 7" to 10". Horizontally and vertically rotatable and tiltable, with theft protection.

Ref. no.

20849

Mount for tablet/iPad



for fixing at rear of a shelf, only mountable ex works. Anodised aluminium, adjustable from 7" to 10". Horizontally and vertically rotatable and tiltable, with theft protection.

Ref. no.

20854

Waste collector "tango"



with practical knee operated opening system and plastic casing, white (RAL 9016). Please order 2 standard profile rails, (ref. no. 26046) for fastening to the cart.

Capacity 8 litres

Dim.:
220 x 220 x 360 mm (w x d x h)

Ref. no.

16701

Capacity 20 litres

Dim.:
350 x 250 x 520 mm (w x d x h)

Ref. no.

16702

Waste collector "tomo"



made of smooth and easy to clean plastic (100% polypropylene), grey (RAL 7045), with detachable, well closing lid and with mount. Capacity of 10 litres (for disposable bags approx. w 450 x h 350 mm).

Dim.:
315 x 180 x 315 mm (w x d x h)

Ref. no.

14097

for **standard profile rail**

Ref. no.

14066

Injection sets / KOMBI-Set / Dressing dispenser / Universal holder /
Reanimation board / Infusion bottle holder / Endoscope holder

**Injection set
"PicBox®"**



for **cart 60**

The dispenser with **first-in/first-out guarantee**, to be fixed between the profiles. Grey-white plastic casing with clear glass boxes and dust protection lids. Consisting of syringe dispenser with 3 boxes and canula dispenser with 5 boxes. To be fixed between the profiles.

Dim.:
600 x 220 x 410 mm (w x d x h)

Ref. no.

26099

KOMBI set



Recommended for taking blood samples and giving injections. Consisting of "TAKEFIX" swab dispenser, "TIPPFIX" swab moisturiser and additional container. All parts are made of unbreakable plastic, grey-white.

Dim.:
120 x 230 x 265 mm (w x d x h)

Ref. no.

13916

**Injection set
"COMPACT"**



for **cart 60**

High quality casing, white plastic veneered, dust-proof sealed by plexiglass covers. Ideally arranged for easy access. To be fixed between the profiles.

Dim.:
600 x 220 x 450 mm (w x d x h)

Ref. no.

13914

**Dressing dispenser
"COMPACT"**



for **cart 60**

High quality casing, white plastic veneered, dust-proof sealed by plexiglass covers. Ideally arranged for easy access. To be fixed between the profiles.

Dim.:
600 x 220 x 450 mm (w x d x h)

Ref. no.

13913

Universal holder "Snap"



made of white plastic, for boxes of disposable gloves, tissues or wipes. The spring-loaded clamp is adjustable to fit many box sizes up to about 1 kg; for single handed operation with mount for doppio.

Ref. no.

26098

**Reanimation board
with holder**



blue, made of sturdy polyethylene, 2 kg lightweight supports up to 150 kg. Holder made of stainless steel. Suitable for trolleys with minimum 380 mm high rear walls.

Dim.:
600 x 420 x 60 mm (w x d x h)

Ref. no.

17370

**Infusion bottle
holder**



for single-handed operation, chrome-plated with 4 hooks, Ø 25 mm.

Ref. no.

26057

**Universal
endoscope holder**



made of plexiglass with several openings, height-adjustable rod for single-handed operation.

Ref. no.

26048

Endoscope tube of plexiglass, open at the bottom, 120 cm long, inside Ø 3,2 cm, outside Ø 4 cm, with holding clips.

Ref. no.

13977

As above, closed at the bottom for disinfectant solution.

Ref. no.

13978

Holder for dispenser boxes / Container for canulas / Ampoule opener / Soap and disinfectant dispenser / Drip tray / Waste container / Disposal boxes

Holder for dispenser boxes



made of brushed stainless steel. Suitable for 1 to 4 common dispenser boxes of disposable gloves or tissues, with mounting material for **doppio**.

Holder 1-fold (for 1 box)
Inner dim.: 142 x 94 x 250 mm (w x d x h)

Ref. no.

13026

Holder 3-fold (for 3 boxes)
Inner dim.: 253 x 94 x 394 mm (w x d x h)

Ref. no.

13027

Holder 4-fold (for 4 boxes)
Inner dim.: 253 x 94 x 525 mm (w x d x h)

Ref. no.

05071

Container for canulas "Schluckfix"



Disposable container for approx. 500 used canulas, made of black recycled plastic.

Dim.:
55 x 100 x 170 mm (w x d x h)
Set with 20 boxes and holder

Ref. no.

13917

Original carton

with 100 boxes, without holder

Ref. no.

09490

Holder

separately for "Schluckfix"

Ref. no.

13918

Ampoule opener "CUPFIX"



with grey-white plastic casing. Receptacle for used ampoules, removable. Diamond file spring-loaded in a stainless steel plate.

Dim.:
70 x 150 x 200 mm (w x d x h)

Ref. no.

13921

Spare diamond file

for "CUPFIX" ampoule opener

Ref. no.

04663

Holder

separately for "CUPFIX"

Ref. no.

13922

Soap and disinfectant dispenser



made of naturally anodised aluminium, for standard bottles of 1 litre, with fixing material.

Dim.:
95 x 185 x 290 mm (w x d x h)

Ref. no.

16087

Drip tray

separately, removable, to clip on the soap and disinfectant dispenser.

Ref. no.

16256

Waste container for instruments



with polypropylene bowl, white, can be sterilized up to 134° C, with mount.

Dim.:
300 x 230 x 130 mm (w x d x h)

with slit lid

Ref. no.


13912

without slit lid

Ref. no.

13960

Disposal boxes "Medi-Müll" à 1.5 ltrs.



made of sunshine yellow polypropylene with HAEBERLE-green lid. Offers 3 practical ways for the disposal of all types of canulas and vacuum containers. Ideally suited for scalpels, pipettes, broken glass, swabs etc.

Height 152 mm, Ø Lid 143 mm, Ø base 108 mm.

Set à 20 boxes and mount.

Ref. no.

13919

Carton à 100 boxes, without mount.

Ref. no.

16653

Mount



separately for "Medi-Müll" box

Ref. no.

13920

Arm support / Catheter holder / Holders for oxygen bottle / Multiple socket outlet / Cable and hose holder / Standard profile rails / Clip-holder / Attachment device

Arm support



extra sturdy, black. Padded with integral upholstered foam, with 2 attachment devices.

Ref. no.

13923

Catheter holder



made of stainless steel, 3 compartments, for hanging on standard profile rail.

Dim.:

160 x 86 x 560 mm (w x d x h)

Ref. no.

13959

Please order standard profile rail (Ref. no. 26046) for lateral mounting.

Holder



for oxygen bottles à 5/10 ltrs. with up to Ø 145 mm; upper holder with fixing tape included. White aluminium (RAL 9006), powder-coated. For mounting on the carrier profile.

Ref. no.

16322

Multiple socket outlet 4-fold, standard



4 grounding sockets, light-grey cable 3 m long, with moulded grounding plug, green 2-pin switch, quality housing, anodised aluminium with 2 grey plastic fastening angles.

Dim.:

269 x 44 x 44 mm (w x d x h)

Ref. no.

21693

Further multiple socket outlets you will find in the separate list of electrical accessories.

Additional upper holder

Ref. no.

16324

Holder



for oxygen bottle of 2 ltrs., white-aluminium (RAL 9006). Upper supporting ring Ø 102 mm, lower supporting ring Ø 84 mm. For mounting to the carrier profile.

Ref. no.

13956

Cable and hose holder



with 2 joints, silvergrey aluminium, including assembly clamp. Reach approx. 950 mm.

Ref. no.

13926

Standard profile rail



made of stainless steel, with 2 fixing points for lateral fastening to the carrier profiles. Length 426 mm.

Ref. no.

26046

Standard profile rail



made of stainless steel, with 2 fixing points for rear mounting on horizontal arm. Length 446 mm

Ref. no.

13953

Clip-holder



for securing cables, made of grey plastic. Fastening in the profile grooves.

Ref. no.

19015

Attachment device



for universal use for all attaching parts to the **doppio** cart. Suitable for round pipes Ø 10 mm and Ø 16 mm.

Ref. no.


13929

All prices are understood exclusive of VAT. As an innovative company we reserve the right to amend prices and designs. As for the rest our General Terms of Business in the valid version apply.






WEEE-Reg. no. DE 20235235



Planning sketch for your **doppio** cart



Chassis / castors / load capacity

	Cart	45 <input type="checkbox"/>	60 <input type="checkbox"/>
	Castors Ø100 mm (load up to 120 kg)	with 2 brakes <input type="checkbox"/>	<input type="checkbox"/>
	45 <input type="checkbox"/>	60 <input type="checkbox"/>	
	Castors Ø125 mm (load up to 120 kg)	with 2 brakes <input type="checkbox"/>	<input type="checkbox"/>
	60 <input type="checkbox"/>	with 2 brakes <input type="checkbox"/>	<input type="checkbox"/>
	Castors Ø150 mm (load up to 120 kg)		
	4 electrically conductive castors		<input type="checkbox"/>

Profile height



standard	variable
566 mm <input type="checkbox"/>	0 - 566 mm <input type="checkbox"/>
950 mm <input type="checkbox"/>	0 - 1206 mm <input type="checkbox"/>
1206 mm <input type="checkbox"/>	0 - 1846 mm <input type="checkbox"/>
566 / 1206 mm <input type="checkbox"/>	

(grids = 32 mm)

For optimal stability, each cart with a profile height up from 950 mm has to be equipped with one horizontal arm, ref. no. 26067 or ref. no. 26068

Colour of decorative strips



- grey
- blue

Space for your ideas:

HAEBERLE

mobile workstations

Our versatile product range



Your medical distributor

HAEBERLE GmbH + Co. KG

Breitwiesenstr. 13, DE-70565 Stuttgart

P.O.Box 800524, DE-70505 Stuttgart

Phone +49 (0) 711/78314-0

Fax +49 (0) 711/78314-77

info@haeberle-med.de

www.haerberle-med.de

Version 10/19