

HAEBERLE

We.Move.Design.®

NEW



Sharps bins
Medi-Müll
Product information

Medi-Müll sharps bins - Safe disposal of medical waste

Correct and safe disposal of medical waste is essential for every medical facility. Syringes and cannulas in particular, which are frequently used in everyday medical practice, require special storage and disposal in order to minimise the risk of infection. These specially developed HAEBERLE disposal containers offer numerous advantages for clinics and medical practices.

HAEBERLE Medi-Müll in detail

- Medi-Müll sharps bins in three different sizes - 0.5 l, 1.5 l or 2.5 l.
- The yellow Medi-Müll sharps bins are conical in shape and equipped with a green, firmly lockable lid and convenient handle.
- This high-quality container in a modern design optimises the disposal of sharps in laboratories, hospitals or medical practices.
- Various helpful mechanisms make the disposal of sharps easier for the medical staff.
- The lid features an innovative interim and final locking mechanism.
- The label design makes it easy to see the maximum filling level.
- Accessories - holder for wall-mounting or at the cart system.
- Medi-Müll sharps bins are certified according to: DIN ISO 23907, UN 3291 and TRBA 250.

The HAEBERLE Medi-Müll sharps bins ensure a safe and hygienic working environment in clinics or medical practices: whether wall-mounted, placed on the work surface or in combination with our HAEBERLE equipment cart systems.



Medi-Müll sharps bin 0.5 l

Product description

- Disposal of sharps such as cannulas, vacutainers, scalpels
- Capacity 0.5 l
- Yellow sharps bin with a green lid and handle
- Puncture resistant
- Material: Polypropylen
- Dimensions: Height 151 mm, Ø 83 mm
- Packaging: Box with 187 sharps bins Medi-Müll, space-saving stacked inside each other
- Holder not included

Ref.-no.

4964

Packaging size

Set with 187 Medi-Müll sharps bins



Medi-Müll sharps bin 1.5 l

Product description

- Disposal of sharps such as cannulas, vacutainers, scalpels
- Capacity 1.5 l
- Yellow sharps bin with a green lid and handle
- Puncture resistant
- Material: Polypropylen
- Dimensions: Height 155 mm, Ø 140 mm
- Packaging: Box with 100 / 20 sharps bins Medi-Müll, space-saving stacked inside each other
- Holder not included

Ref.-no.

4965

Packaging size

Set with 100 Medi-Müll sharps bins

4967

Set with 20 Medi-Müll sharps bins



Medi-Müll sharps bin 2.5 l

Product description

- Disposal of sharps such as cannulas, vacutainers, scalpels
- Capacity 2.5 l
- Yellow sharps bin with a green lid and handle
- Puncture resistant
- Material: Polypropylen
- Dimensions: Height 240 mm, Ø 140 mm
- Packaging: Box with 55 sharps bins Medi-Müll, space-saving stacked inside each other
- Holder not included

Ref.-no.

4966

Packaging size

Set with 55 Medi-Müll sharps bins



Medi-Müll sharps bin holder

Product description

- Bracket for mounting on wall or at cart system
- For attaching 1.5 l / 2.5 l Medi-Müll sharps bins
- White aluminium (RAL 9006)
- Powder-coated
- Mounting material included



Ref.-No.	Packaging unit
4671	1 pc. / Wall mounting
4670	1 pc. / Variocar, Variocar Viva and 08/16
13920	1 pc. / swingo and doppio
20232	1 pc. / toro
18059	1 pc. / keo

Medi-Müll holder standard rail

Product description

- For mounting on the left or right of all standard profile rails
- Mounting material included
- Order standard rail separately



Ref.-No.	Packaging unit
16677	1 pc.

Get in touch!

Benefit from the advantages of our new HAEBERLE disposal boxes for a safe and hygienic working environment in healthcare facilities. As a specialist we offer you an economically attractive solution for your customers with the new Medi-Müll sharps bins. HAEBERLE Medi-Müll minimises the risks involved in disposing of medical waste and ensures a high level of hygiene. Find out more and contact us.

Send us an e-mail to service@haeberle-med.de. We look forward to receiving your enquiry.